



# GLOBAL EMERGENCY MEDICAL REGISTRY

## Obtaining Weight using Bed Scale Skill Documentation Form

Candidate (Print): \_\_\_\_\_ Date: \_\_\_\_\_

Examiner (Instructor or Licensed Provider): \_\_\_\_\_

Examiner Signature: \_\_\_\_\_

Pass \_\_\_\_\_ Fail \_\_\_\_\_

Task	Correct	Incorrect
Dons appropriate PPE and explains procedure, speaking clearly, slowly, and directly, maintaining face-to-face contact whenever possible		
Position scale at bedside		
Turn on the scale and roll patient/bed onto scale		
Use the "Tare" or "Preset Tare" function to remove the weight of the bed		
Roll the patient onto the scale		
Read the patient's weight		
Remove patient from scale		
Wash hands after removing PPE		
Records weight based on indicator on scale		

### Critical Failure Criteria

Failure to complete within 10 minute time limit
Failure to take appropriate PPE precautions
Contaminates equipment or site without appropriately correcting the situation
Failure to manage the patient as a competent provider
Exhibits unacceptable affect with patient or other personnel
Uses or orders a dangerous or inappropriate intervention

NOTE: You must factually document any "incorrect" or critical failure criteria on back of this form